

Excellence in Design

Based on the award winning design of the Elite Series Plethysmograph, the new *Platinum Elite* offers comfort, accessibility and room to spare to meet the pulmonary testing needs of all patients, pediatric through adult.



Ideal for Small Testing Rooms and Crowded Labs

- Small 42 inch diameter foot-print
- Unique "zero clearance" door

Easy Access for all Patients

- Patients enter the cabin without the need to "step-up" and are seated with minimal effort
- Patients can be transferred from a wheelchair with ease, or testing can be performed outside the cabin



Relaxed and Pleasant Patient Testing Environment

- Direct unobstructed patient view
- Expanded interior space comfortably accommodates large and tall patients without patient weight limitations
- Inflatable door seal with no mechanical latches provides additional safety and comfort